



Securely email genetic reports via Sonic Dx




Sonic Dx uses best-practice encryption to safely transfer files or share reports with your patient or other healthcare providers, to facilitate patient management and follow-up.


- 1 Open Sonic Dx on your desktop or handheld device. Click on 'Search' and enter the laboratory request number or the patient's name. Please ensure that the 'Last X Days' field is adjusted to include the date of the test.




- 2 Open the report.



- 3 Click on the 'Share' option on the right hand side of the 'genetic report' heading, e.g. NIPT, PGx, FH etc.



- 4 Enter the recipient's email address, a PIN and their mobile phone number (optional). The recipient will need this PIN to access their copy of the report.



- 5 Sonic Dx will send the recipient a link to the PDF report, which they can securely download after entering the PIN created in Step 4.

Please note this guide relates to Sonic Dx as it appears on desktop and handheld devices will vary.