

Beacon expanded carrier screening

Practice management software template for Genie


How to request Beacon expanded carrier screening through Genie

We are pleased to be able to provide a built-in template for our Beacon expanded carrier screening request form.

This means you can now access the Beacon expanded carrier screening form within Genie, populate with patient details and print off a request form, as you would with standard pathology requests.


Step 1

Open patient file.




Step 2

Click on letter icon in the patient file.



Step 3

Select expanded carrier screening and open (please note if template is not stored, please follow adjacent instructions).



To download template:

Download the template first via Special
→ Software Updates → Forms →
Request → Beacon Expanded Carrier
Screening Test (Sonic).

This is only required to be performed once.

Correct at time of printing. Any changes to the software made by the vendor may result in changes to the above procedure for utilising this template. For further assistance, please contact Client IT at your local Sonic Healthcare pathology provider.